

# Photron FASTCAM Analysis (PFA)

## Ver.1.1 Release Notes

Ver.1.1.0.1

Release Date: 19<sup>h</sup> August, 2013

### Changes

---

- ✓ License agreement is modified.

Ver.1.1.0

Release Date: 23<sup>rd</sup> July, 2013

### New features

---

- ✓ Magnified images of each tracker.
- ✓ Export a movie file on which a tracked image and a graph image are displayed.
- ✓ Export a project file to resume PFA.
- ✓ On- line activation.
- ✓ Multiple language support. (English, Japanese, Chinese)
- ✓ 64-bit installation.
- ✓ Plots on sample points on graphs
- ✓ Zoom control for graphs to change the position and zoom rate.
- ✓ Control of current frame on the graph.
- ✓ Load an image file which is not opened from a cih file in PFV.
- ✓ Display a cross cursor when user can add a tracker.
- ✓ Selection of trackers by dragging an area and deletion of selected trackers from "Delete" key on keyboard are enabled.
- ✓ Change of trackers position by dragging directly on a failure frame is enabled.

### Changes

---

- ✓ User interfaces are improved.

- ✓ The icon is changed to the original one.
- ✓ Windows 8 is supported.
- ✓ The license system is changed to a newer version.
- ✓ Trial period is changed to 90 days.