

# Photron FASTCAM Analysis (PFA)

## Ver.1.2 Release Notes

Ver.1.2.0

Release Date: 16<sup>th</sup> June, 2014

### New features

---

- ✓ Standalone launch.
- ✓ New tracking algorithm.(Binarized center of gravity)
- ✓ Tracking point can be adjusted at the zoom window.
- ✓ Scale is added on over view graph area.
- ✓ Jump to frame from a graph point.
- ✓ Selectable output files.

### Changes

---

- ✓ Windows 8.1 is supported.
- ✓ Bayer files which are generated from FASTCAM Mini UX50 become available.
- ✓ MRAW, RAW, RAWW files which are generated from FASTCAM Mini UX50 become available.
- ✓ Tracking point is changed to a cross mark from a circle mark.
- ✓ Zoom window is fitted to template size.
- ✓ Zoom window refresh during tracking is revised to smooth.
- ✓ Time unit on the graph is changed from “ms” to “s”.
- ✓ No license PC's functional limitation is changed.
- ✓ Language list is revised.